

| | | | | | |
|-------------------|---|-----------|--------------|-------|-----------------|
| primotec ® | Material Safety Data Sheet Separate Product / Preparation | SDB – No. | 04-01 | Page | 1 of 2 |
| | | Revision | 1.0 | Date: | 7/1/2004 |

| | | | | | | | |
|----------|------------------------------------|-----------------------|--|-----------|--------------------|-----------------------|------------------------|
| 1 | Material and corporate name | | | | | | |
| 1.1 | Details on the product | Trade Name | Metaseal separator | | | Item No. IM300 | |
| 1.2 | Commercial enterprise | Name | primotec, Joachim Mosch | | | | |
| | | Address | Tannenwaldallee 4, D-61348 Bad Homburg | | | | |
| | | Telephone | 0049-6172-99770-0 | Fax | 0049-6172-99770-99 | e-mail | primotec@primogroup.de |
| | | Information about | Safety | Telephone | 0049-6172-99770-0 | Fax | 0049-6172-99770-99 |
| | | Emergency information | See above | | | | |

| | | |
|----------|---|-----------------------------|
| 2 | Composition / information on ingredients | |
| 2.1 | Description | Emulsion, pigment and water |

| | | |
|----------|---|----------------|
| 3 | Hazards identification | |
| 3.1 | Special hazards information for human and environment | Not applicable |

| | | |
|----------|---------------------------|---|
| 4 | First aid measures | |
| 4.1 | General Information | No special measures necessary |
| 4.2 | After contact with skin | Wash hands thoroughly with water and soap, consult a medical doctor in case of irritation |
| 4.3 | After contact with eyes | Rinse thoroughly with much water, seek medical advice, if irritation |
| 4.4 | After ingestion | Cause vomiting and consult a medical doctor |

| | | |
|----------|---|--|
| 5 | Fire fighting measures | |
| 5.1 | Suitable extinguishing media | |
| | Powder water spray extinguisher | |
| 5.2 | Unsuitable extinguishing media for safety reasons | |
| | Water jet | |
| 5.3 | Special productive equipment for fire fighting | |
| | Not applicable | |

| | | |
|----------|--|--|
| 6 | Accidental release measures | |
| 6.1 | Precautionary measures related to people | |
| | Not applicable | |
| 6.2 | Environmental protective measures | |
| | Not required | |
| 6.3 | Methods of cleaning / adsorption | |
| | Absorb with cloth or similar tissue and throw away | |

| | | |
|----------|---|--|
| 7 | Handling / Storage | |
| 7.1 | Instructions on safe handling | |
| | No special measures necessary | |
| 7.2 | Information on fire and explosion protection | |
| | No special measures necessary | |
| 7.3 | Storage | |
| | Keep at a cool and dark place. Not above 20°C | |
| 7.4 | Further information about storage | |
| | Not required | |

| | | | | |
|----------|--|---|-----------------|---|
| 8 | Exposure controls / personal protection | | | |
| 8.1 | Additional information on shaping of technical constructions | | | |
| | No further information | | | |
| 8.2 | Exposure limits | | | |
| | None | | | |
| 8.3 | Personal protection | | | |
| | No specific recommendations | | | |
| | Breath protection | Not required | Hand protection | Allergic persons are required to wear non rubber gloves |
| | Eye protection | Don't let the product come into contact with the eyes | Body protection | Wear common industrial safety clothing |
| | Hygiene measures | Keep common hygiene for the job | | |

| | | | | | |
|-------------------|---|-----------|--------------|-------|-----------------|
| primotec ® | Material Safety Data Sheet Separate Product / Preparation | SDB – No. | 04-01 | Page | 2 of 2 |
| | | Revision | 1.0 | Date: | 7/1/2004 |

| | | | | | | | |
|----------|---|------|----------------|--------------------------|-------|----------------|---------|
| 9 | Physical and chemical properties | | | | | | |
| 9.1 | Appearance: | Form | Liquid | Color | Clear | Odor | Typical |
| 9.2 | Specifications pertinent to safety | | | | | | |
| | Burn out temperature | | Not determined | Burn out ash measurement | | None | |
| | Flashpoint | | Not determined | Self inflammation | | Not applicable | |
| | Danger of explosion | | None | Explosion limit | | Not applicable | |
| | Vapor pressure | | Not determined | Density | | Not determined | |
| | Solubility in water | | Soluble | Ph-value | | Not applicable | |

| | | |
|-----------|----------------------------------|--|
| 10 | Stability and reactivity | |
| 10.1 | Thermal decomposition | No decomposition when used as directed |
| 10.2 | Hazardous reaction | None |
| 10.3 | Hazardous decomposition products | None when used as prescribed |

| | | |
|-----------|----------------------------------|--------------|
| 11 | Toxicological information | |
| 11.1 | Skin irritation | Non-irritant |
| 11.2 | The primary irritation index | None |

| | | |
|---|-------------------------------|--|
| 12 | Ecological information | |
| No information about negative ecological consequences because of the products consistency | | |

| | | |
|-----------|---|---|
| 13 | Disposal consideration (Germany) | |
| 13.1 | Product | Observe local/national decrees; usually no special treatment; don't give into water |
| 13.2 | Package | The package is made of plastic and is disposed as household waste; the box can be reused after cleaning |

| | | |
|---|--|--|
| 14 | Packing and transport decrees (Germany) | |
| No special packing and transportation necessary because of non hazardous material | | |

| | | |
|-----------|--|----------------------------|
| 15 | Miscellaneous decrees (Germany) | |
| 15.1 | Licence according to GefStoffV/EG | No duty |
| 15.2 | National decrees | None |
| 15.3 | Grade of endangering water | Grade 1 (self- assessment) |

NOTE

This Safety Data Sheet is based on our current knowledge of and experience with the product. It is exclusively meant to reduce risk factors handling the product. Hence it only describes aspects of safety. The details are no guaranteed characteristics according to our public terms of business. They do not release the user from the compliance of laws, regulations and rules.

Our Safety Data Sheets are free of charge with the first subscription. Additional copies can be requested for refund of expenses. Updating the Safety Data Sheets in accordance with the relevant regulations, we inform all users, who received the Safety Data Sheet in the last 12 months before the date of revision. This occurs by sending the revised version free of charge. Owners of out-dated Safety Data Sheets can enquire about the particular level of revision.